

**Product data sheet**  
**Mounting plate with cover 950.50.10040**

# HEWI




**Mounting plate with cover 950.50.10040**

- rectangular mounting plate with matching, rectangular cover
- used to hold the mobile hinged support rails
- 130 mm wide, 181 mm high, 15 mm deep
- mounting plate made of stainless steel
- for mounting with non-corrosive and tested HEWI fixing materials for different wall systems (to be ordered separately)
- cover frame made of high-quality polyamide in HEWI colour 98 (signal white)
- face plate made of mirror polished stainless steel
- CE marking: Medical devices Class 1 according to Directive 93/42 / EEC



Dimensions in mm



Surface/Colour

XC   mirror polished	98
Stainless steel	Polyamide

Technical information subject to alteration, 15.06.2015

HEWI Heinrich Wilke GmbH  
Postfach 1260  
D-34442 Bad Arolsen

Telefon: +49 5691 82-366  
Telefax: +49 5691 82-99366  
international@hewi.com  
www.hewi.com

HEWI (UK) Limited  
Holm Oak Barn, Beluncle Halt  
Stoke Road  
Hoo Kent ME3 9NT

Phone +44 1634 258 200  
Fax +44 1634 250 099  
info@hewi.co.uk  
www.hewi.co.uk