


Product data sheet

Hinged support rail Mono 950.50.61040

# HEWI




### Hinged support rail Mono 950.50.61040

- with a horizontal, two-dimensional gripping level
- for holding on to and providing support
- rail made of rectangular profile, 40 mm wide and 30 mm, provides a wide support surface and large radii for pleasant gripping
- 600 mm projection, wall mounting plate 181 mm high and 130 mm wide
- can be folded upwards and, folded downwards with braking mechanism
- made of metal, high-quality chrome-plated
- wall-mounting plate cover made of high-quality chrome-plated synthetic material
- stable three-point fixing
- for wall mounting with concealed, non-corrosive and tested fixing materials for different wall systems
- with toilet roll holder 950.50.01150, toilet roll holder with WC flushing mechanism and function button (radio-controlled) 950.50.06490 (right), 950.50.06590 (left) and arm pad 950.50.03192 for retrofitting
- CE marking: Medical devices Class 1 according to Directive 93/42 / EEC
- fulfils the DIN 18040 requirements



#### Dimensions in mm



Surface

chrome

Technical information subject to alteration, 15.06.2015

HEWI Heinrich Wilke GmbH  
Postfach 1260  
D-34442 Bad Arolsen

Telefon: +49 5691 82-366  
Telefax: +49 5691 82-99366  
international@hewi.com  
www.hewi.com

HEWI (UK) Limited  
Holm Oak Barn, Beluncle Halt  
Stoke Road  
Hoo Kent ME3 9NT

Phone +44 1634 258 200  
Fax +44 1634 250 099  
info@hewi.co.uk  
www.hewi.co.uk