



**Product data sheet**

**Hinged seat Premium White Edition 802.51.W220**

**HEWI**




**Hinged seat Premium White Edition 802.51.W220**

- ergonomically shaped seat with back and armrests
- used for safe sitting in the shower area
- seat weight capacity up to 150 kg
- seat and backrests easy to fold up and down (lock in upper position)
- space saving
- 620 mm wide, 420 mm high, 455 mm deep
- seat 444 mm wide, 380 mm deep, tilt adjustable
- seat height adjustable to four levels from 420 mm to 510 mm (every 30 mm)
- easy to install due to separate fixing of HEWI wall plate and seat
- concealed, non-corrosive HEWI fixing material for different wall systems
- seat and backrest made of high-quality PUR in the two HEWI colours 98 (signal white) and 99 (pure white)
- seat with hygienic opening and water run-off
- hinge made of white powder-coated aluminium with highly polished chrome design features
- armrests made of high-quality PUR in the two HEWI colours 98 (signal white) and 99 (pure white)
- armrests can be folded up and down separately
- CE marking: Medical devices Class 1 according to Directive 93/42 / EEC
- fulfils the DIN 18040, ÖNORM B1600/1601 and SIA 500 requirements



**Dimensions in mm**



**Surfaces/Colour**

chrome	white	98	99
--------	-------	----	----

**Alternative surfaces/colours**

chrome	matt	97	92	98	99	97	92
--------	------	----	----	----	----	----	----

Arm rests

Technical information subject to alteration, 15.06.2015

HEWI Heinrich Wilke GmbH  
Postfach 1260  
D-34442 Bad Arolsen

Telefon: +49 5691 82-366  
Telefax: +49 5691 82-99366  
international@hewi.com  
www.hewi.com

HEWI (UK) Limited  
Holm Oak Barn, Beluncle Halt  
Stoke Road  
Hoo Kent ME3 9NT

Phone +44 1634 258 200  
Fax +44 1634 250 099  
info@hewi.co.uk  
www.hewi.co.uk